

INSTRUCTION AND INFORMATION ABOUT THE MANUFACTURER

nTr Professional Faceshield - consists of extremely lightweight and elegant frame in which through a click system a reusable screen is installed. The frame is placed on the head with the help of an eraser that can be adjusted individually. For additional comfort a mellow protector for the forehead is placed. The set consists of:

Frame - 1p. made of PP with protector of EVA and eraser with width 10.5mm

Front Protective Window – 5p. that are 0.2mm thick and made of PET Instruction and Information for the Producer – 1p. Box - 1p.

nTr Professional Faceshield – protects from drops and splashes of liquids, responding to the European Regulation (EC) 2016/425, European Standard EN166:2001 point 7.2.4 b, EN168:2008 objection 12.2 and the requirements for EN166:2001 point 6.1, 6.2, 6.3;

nTr Professional Faceshield – can be used with mask on as well as with protective or optical glasses.

nTr Professional Faceshield – is incredibly lightweight – only 46 g. it does not disturb the eye vision or the free movement of the head. It can be worn for hours without fatigue.

nTr Professional Faceshield – is not disinfected. Before usage disinfect the faceshield. Suitable for disinfection are all disinfectors containing ethanol. Avoid solvents, ketones, strong acids and alkaline elements.

nTr Professional Faceshield – store in a dry place with temperature from -25C to +50C in its box for no longer than 5 years.

Manufacturer : TRANSGRAPH LTD, 1528 Sofia, BULGARIA, str. Nedelcho Bonchev 10, www.transgraph.eu , shield@transgrapf.eu